BeneCheck™ Uric Acid

BK-U1

For use with the BeneCheck Meter PD-G001-3-P/BSU010/BSU01D.

IMPORTANT:
Please read this instruction and user's manual before using BeneCheck Uric Acid Test Strip (Model: BK-U1).

■ Intended Use
- Measuring uric acid level in capillary whole blood from fingertip.
- The meter can be used by laypersons or healthcare professionals.
- *In vitro* diagnostic use only.

■ Principle of the Test
- The blood sample is drawn into the strip reaction zone by capillary force. The electrodes oxidize uric acid in whole blood and activate a current that is proportional to uric acid in blood. Then the current is converted into a reading of uric acid concentration.

■ Chemical Composition of Strip
- Wetting agent: cellulose
- Other ingredients (buffer, surfactant, etc.): ≥ 0.01 mg

■ Storage and Handling
- Storage the strips within 4-30 °C (39-86 °F), use within 3 months after first open.
- Handle the strips with dry and clean hand.
- Recap strip vial immediately after use.
- Use the strip immediately after taking out from the vial.
- Do not use expired strips. It may cause false result.
- Do not expose the strips under sunlight, heat or humid environment.
- Do not repack, change vial, refrigerate or freeze the strips.
- Do not bend, cut or fold the strips.

■ Auto Strip Recognition
Please refer to User's Manual for more information.

■ Coding Procedure
Insert the code strip → Code number and test mode display on screen → Make sure the code on screen/ code strip/ strip vial label is the same → Make sure test mode on top of screen display "UA"

■ How to Perform a Test
Take a strip → Recap strip vial immediately → Insert the strip → Make sure the code is correct → Sampling by lancing device → Touch the blood sample with strip → Result will show after countdown → Discard used strip and lancet to waste container
Please refer to User's Manual for detail procedure.

IMPORTANT: Please double check the test mode before testing.

■ Test Results
- Measuring range: 3-20 mg/dL
- The normal range for healthy male: 3.4-7.0 mg/dL
 The normal range for healthy female: 2.4-6.0 mg/dL
 (The range is only a reference.)
- If "Lo" is displayed on the meter, uric acid level may be below 3 mg/dL.
- If "Hi" is displayed on the meter, uric acid level may be over 20 mg/dL.

If the test result is unusually low, high or inconsistent with how you feel, please check the follows and repeat the test again with a new strip and make sure you have followed the instructions carefully.
- Strips are not expired.
- Meter has been coded properly, and the code number is consistent with the strip.
- Meter and strips are working properly.
- Blood sample has filled up the strip inlet sufficiently.

If you still have doubt with the test result, please consult healthcare professional or contact authorized distributor.

IMPORTANT:
- Please consult healthcare professional before making medical decision based on the tested results.
- Strips from other brand are not compatible to BeneCheck's meter.
- POTENTIAL BIOHAZARD. The used test materials could be sources of infection, please beware especially when the system is used by multiple users.

■ Control Solution Test
Control solution is used to check the performance of the system. If the control solution test result falls within the range listed on the strip vial label, the system is working properly. Control solution is not included in the system. Please contact an authorized distributor for purchasing.
A control test should be performed:
- When the system does not work properly.
- When the test result is unusual or inconsistent.

Out of range control solution test result may be caused by:
- Performing the test inaccurately.
- Contaminated or expired control solution.
- System deterioration or malfunction.
- The codes between meter and strip are inconsistent.

IMPORTANT: Each strip vial might have different control range. Please check every time when you perform a control solution test. Do not use the system if the control solution tests keep result wrongly. Please contact authorized distributor for assistance.

■ Limitations
For getting an accurate test result, please read the following limitation carefully.
- **Single use only. Do not reuse the strips.**
- Do not perform the test under strong electronic field with signal interference.
- Do not perform the test on newborns.
- The test strip can be used at altitude up to 10,000 feet (3048 meters).
- Hematocrit (Hct) should be between 30-55 %.
- Use only fresh capillary whole blood from fingertip.
- Do not use serum, plasma or venous whole blood.

■ Additional Information for Healthcare Professionals
- When using venous whole blood for tests, please collect with heparin-containing test tube within 30 mins after drawing. Test tubes contain EDTA or Sodium Fluoride may affect test result.
 Test tubes contain Sodium Citrate will not affect test result.
- The 23 interference substances result in insignificant bias to the measurement of uric acid concentration within the therapeutic range. Refer to CLSI EP07 and EP37 for additional reference source of the therapeutic range.

Acetaminophen	3 mg/dL	Glucose	503 mg/dL
Ascorbic acid	1.5 mg/dL	Ibuprofen	40 mg/dL
Allopurinol	10 mg/dL	Indomethacin	20 mg/dL
Amiloride	1.0 mg/dL	Ketoprofen	20 mg/dL
Atenolol	10 mg/dL	Metformin	40 mg/dL
Bilirubin	5 mg/dL	Methyldopa	1.0 mg/dL
Creatinine	10 mg/dL	Salicylate (Aspirin)	50 mg/dL
Cholesterol	322 mg/dL	Tetracycline	6 mg/dL
Colchicine	10 mg/dL	Tolazamide	25 mg/dL
Diclofenac	5 mg/dL	Tolbutamide	40 mg/dL
Dopamine	1.6 mg/dL	Triglycerides	1500 mg/dL
Glibenclamide	80 mg/dL		

- Gentisic acid at therapeutic levels (> 2.5 mg/dL) may affect the test result.
- Decreased peripheral blood flow may case false results. Ex: severe dehydration, in shock, or in a hyperosmolar state (with or without ketosis), hypotension.
- The test strips are not recommended for critically ill patients.

■ Accuracy:
The accuracy of the uric acid test was performed with capillary finger blood and compared with reference instrument calibrators, which were referenced to the NIST SRM 913 Standard.

At uric acid concentration < 5 mg/dL

	Within ± 0.25 mg/dL	Within ± 0.5 mg/dL	Within ± 0.75 mg/dL	Within ± 1 mg/dL
PD-G001-3-P	49/66 (74.2 %)	54/66 (81.8 %)	65/66 (98.5 %)	66/66 (100.0 %)
BSU010, BSU01D	51/90 (56.7 %)	78/90 (86.7 %)	90/90 (100 %)	-

At uric acid concentration ≧ 5 mg/dL

	Within ± 5 %	Within ± 10 %	Within ± 15 %	Within ± 20 %
PD-G001-3-P	210/234 (89.7 %)	214/234 (91.5 %)	225/234 (96.2 %)	234/234 (100.0 %)
BSU010, BSU01D	83/210 (39.5 %)	139/210 (66.2 %)	203/210 (96.7 %)	-

■ Precision:
Repeatability:

PD-G001-3-P					
Mean (mg/dL)	3.2	6.9	9.7	12.5	15.7
Mean (mmol/L)	0.19	0.41	0.58	0.74	0.93
SD (mg/dL)	0.2	0.3	0.6	0.6	0.9
SD (mmol/L)	0.01	0.02	0.03	0.04	0.06
CV (%)	-	4.6	5.8	4.8	6.0

BSU010, BSU01D					
Mean (mg/dL)	4.1	6.2	9.9	12.9	18.3
Mean (mmol/L)	0.24	0.37	0.59	0.77	1.09
CV (%)	5.2	4.7	4.6	4.2	3.8

Intermediate precision:

PD-G001-3-P			
Mean (mg/dL)	4.3	7.4	12.0
Mean (mmol/L)	0.25	0.44	0.72
SD (mg/dL)	0.3	0.4	0.6
SD (mmol/L)	0.02	0.03	0.04
CV (%)	-	5.7	5.1

BSU010, BSU01D			
Mean (mg/dL)	4.60	7.20	13.10
Mean (mmol/L)	0.273	0.428	0.779
SD (mg/dL)	0.16	0.21	0.31
SD (mmol/L)	0.010	0.013	0.018
CV (%)	3.6	3.0	2.4

■ Symbols
Please refer to User's Manual for more information.

■ Manufacturer
General Life Biotechnology Co., Ltd.
6, 7F., No. 669, Jhongjheng Rd., Shin Juang Dist., New Taipei City 242, Taiwan
Tel：+ 886-2-22053698
Website：www.BeneCheck.com.tw
Email：service@glbiotech.com.tw

Importer/Service/Distribution in Germany:
PHARMADOC GmbH
Hauptstr. 27F, 23923 Lüdersdorf
TEL: 038821-62080
FAX: 038821-62078
E-Mail: mail@pharmadoc.info

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Pharmadoc GmbH
EC REP MDSS GmbH