beurer EM 10 Mini-Pad BODY

CE 0483

Beurer GmbH,
Söflinger Str. 218,
89077 Ulm, Germany
www.beurer.de
Mail: kd@beurer.de

DEUTSCH

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ENGLISH

Please read these instructions for use carefully and keep them to refer to later. Be sure to make them accessible to other users and observe the information they contain.

Please read these instructions carefully before using the device.

LOT	LOT number
	Expiry date
	Manufacturer
3V Typ CR2032	1 x CR2032 battery
	Important! Read the instructions for use.
CE	This device complies with EU Directive 93/42/EEC

1. Getting to know your device

1.1 What is the mini pad?

Migraines, cramps, acute pain, pain resulting from injury and chronic pain are all natural warning signs and need immediate relief. The innovative mini pad is a self-adhesive TENS (transcutaneous electrical nerve stimulation) device for pain relief. The Mini pad uses TENS/EMS technology for targeted pain relief. Made from medical-quality silicone rubber, the elegantly designed Mini pad is extremely slim and flexible and adapts perfectly to the area to be treated. The high-tech circuit provides energy-efficient treatments of 20 minutes each (up to 20 treatments per battery).

Mini pad is a safe and easy way of relieving pain and offers users a completely new experience in dealing with pain.

1.2 What is TENS and what can it do?

TENS, **or transcutaneous electrical nerve stimulation**, relates to the electrical stimulation of the nerves via the skin. TENS is an effective non-pharmacological method of treating different types of pain that have a variety of causes and has no side-effects if administered correctly. The method has been clinically tested and approved and can be used for simple self-treatment. The pain-relieving or pain-suppressing effect is achieved by inhibiting the transference of pain to nerve fibres (caused mainly by high-frequency impulses) and by increasing the secretion of endorphins in the body. Their effect on the central nervous system reduces the sensation of pain. The method is scientifically substantiated and approved as a form of medical treatment. Any symptoms that could be relieved using TENS must be checked by your GP. Your GP will also give you instructions on how to carry out a TENS self-treatment regime.

TENS is clinically tested and approved to treat the following complaints:
- Back pain, particularly in the lumbar/cervical spine area
- Sore joints (e.g. knee, hip and shoulder joints)
- Neuralgia
- Menstrual cramps in women
- Pain resulting from injury to the musculoskeletal system
- Pain caused by circulatory disorders
- Chronic pain with various causes

TENS devices work by passing electrical currents over the skin. The gel film is used as a transfer medium and is subject to natural wear and tear. The gel film must be replaced if it stops providing sufficient contact, as this will prevent the mini pad from sticking to the skin. If it is not replaced, the partially increased current density could irritate the skin.

1.3 Included in the delivery
- 1 x mini pad
- 1 x self-adhesive gel film
- 1 x 3V CR2032 battery

Replacement items
- 1 x self-adhesive gel film and 1 x battery, item no.: 647.15

When necessary, please contact customer services at the address provided in good time.

IMPORTANT!
The mini pad may only be used as follows:
- On people.
- For the intended purpose and as specified in these instructions for use. Improper use can be dangerous.
- For external use only.
- Use only with the original accessories supplied, which can be re-ordered. Failure to do so invalidates the warranty.

2. Important information

The device is not a substitute for medical consultation and treatment. Consult your doctor first if you are experiencing any pain or are suffering from an illness.

WARNING!
To avoid damage to health, we strongly advise against using the mini pad in the following situations:
- If you have implanted electrical devices (e.g. a pacemaker).
- If you have metal implants.
- If you use an insulin pump.
- If you have a high fever (e.g. > 39 °C).
- If you have a known or acute cardiac arrhythmia, or disorders of the heart's impulse and conduction system.
- If you suffer from a seizure disorder (e.g. epilepsy).
- If you are pregnant.
- If you have cancer.
- After an operation, if strong muscle contractions could affect the healing process.
- The device must never be used near to the heart. The stimulation electrodes must not be placed on any part of the front ribcage (where the ribs and breastbone are located), especially not on the two large pectorals. This can increase the risk of ventricular fibrillation and induce cardiac arrest.
- On the skeletal skull structure, or around the mouth, throat or larynx.
- In the neck area/carotid artery area/above the neck.
- In the genital area.
- On acutely or chronically diseased (injured or irritated) skin (e.g. inflamed skin – whether painful or not, reddened skin, rashes, e.g. allergies, burns, bruises, swellings, both open and healing wounds, and post-operative scars where the healing process could be affected).
- In humid environments (e.g. in the bathroom) or when bathing or showering.
- Do not use after consuming alcohol.
- Do not use if connected to a high-frequency surgical device.

Before using the device, consult your doctor if any of the following applies to you:
- Serious illnesses, in particular if you suspect or have been diagnosed with high blood pressure, a blood coagulation disorder, propensity to thrombo-embolic conditions or recurrent malignant growths.
- Skin diseases or open wounds.
- Unexplained chronic pain in any part of the body.
- Diabetes.
- Any sensory impairment that reduces the feeling of pain (e.g. metabolic disorders).
- If you are receiving medical treatment.
- If you are already having stimulation treatment for existing conditions.
- If you suffer from persistently irritated skin due to long-term stimulation at the same electrode site.

PRECAUTIONS:
- Hold the device away from sources of heat and do not use it in close proximity (c. 1 m) to shortwave or microwave devices (e.g. mobile phones), as doing so can result in unpleasant current peaks.
- Do not expose the device to direct sunlight or high temperatures.
- Protect the device from dust, dirt and humidity. Never immerse the device in water or other liquids.
- The device is suitable for self-treatment.
- If the device does not work properly, or you feel unwell or experience pain, stop using it immediately.
- Switch off the mini pad before you attach or remove it.
- Do not adapt the mini pad in any way, for example, by cutting or trimming it. This will produce a higher current density and can be dangerous.
- Do not use whilst asleep, driving a vehicle or operating machinery.
- Do not use whilst undertaking any activity where an unexpected reaction (e.g. strong muscle contractions even at low intensity) could be dangerous.
- Ensure that no metallic objects (e.g. belt buckles or necklaces) come into contact with the electrodes during stimulation. If you are wearing jewellery or have piercings in the area to be treated (e.g. a navel piercing), these must be removed before using the device. Failure to do so could result in spot burns.
- Keep the device away from children.
- Pre-treating the skin with greasy creams or lotions is not recommended.

Damage
- If the mini pad is damaged or you have any doubts, do not use it and contact your retailer or customer services at the address provided.
- Switch the device off immediately if it is faulty or not working properly.
- Do not attempt to open and/or repair the mini pad yourself. Repairs may only be carried out by customer services

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