FORA® 6

ver 2.0 2024/04
312-4634100-019

β-Ketone Test Strip

EN

Warnings

► For *in vitro* diagnostic use (for use outside of the body only). ► For single use only. ► Healthcare professionals and other users testing multiple patients with this systemshould handle everything that has comes into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects. ► Please read this sheet and your FORA 6 Duo / FORA 6 Plus / FORA 6 Connect Multi-functional meter Owner's Manual before you use this test strip. ► Use only FORA 6 β-Ketone test strips with FORA 6 Duo / FORA 6 Plus / FORA 6 Connect Multi-functional meter to obtain accurate results, and be covered by the manufacturer's warranty. ► Results may be inaccurate when testing on patients with abnormally low bloodpressure or those who are in shock. ► Please do not use FORA 6 Duo / FORA 6 Plus / FORA 6 Connect Multi-functional monitoring system on critically ill patients. The collection of capillary blood from the approved sample sites is not advised when the peripheral circulation is impaired as the delivery of physiological β-Ketone level might not be a true reflection. The following circumstances may apply: severe dehydration attributable as a result of diabetic ketoacidosis or due to stress hyperglycemic, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease. ► Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

Intended Use

FORA 6 β-Ketone test strips, when used together with FORA 6 Duo / FORA 6 Plus / FORA 6 Connect Multi-functional meter, allow your β-Ketone levels to be measured by yourself at home or by healthcare professionals. It uses fresh whole blood samples from the finger. This system is not intended for use in the diagnosis or screening of diabetes mellitus.
Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing. Use ONLY heparin for anticoagulation of whole blood samples. Please do NOT use EDTA for anticoagulation.

Contents

► β-Ketone Test Strip ► Code strip ► Instruction for use

Limitations

► Hematocrit: The hematocrit level is limited to between 10% and 70%. Please ask your healthcare professional if you do not know your hematocrit level. ► Neonatal Use: This test strip must not be used for the testing of newborns. ► This test strip is used for testing fresh capillary and venous. ► Altitude Effects: Altitudes up to 10,742 feet (3,275 m) do not affect test results. ► Please see Appendix: Summary of substances and concentrations in excess of limitation with interference.


Storage and Handling

IMPORTANT: Do not use the test strips if they have expired.

► Test strips expire 6 months after first opening. Write the first opening date on the test strip vial when you first opened it. (For vial strip only). ► Use each test strip immediately after taking it out of the vial or the individual foil packet. ► Close the vial immediately after taking out a strip. Keep the vial closed at all times. (For vial strip only) ► Store the test strips in their original vial ONLY. Do not transfer them to a new vial or any other containers. (For vial strip only) ► Do not touch the test strips with wet hands. ► Do not bend, cut, or alter the test strip. ► Store the test strips in a cool, dry place between 2°C and 30°C (35.6°F and 86°F) and 10% to 85% relative humidity. ► Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.

Testing Your β-Ketone

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

Reading Your Result

The β-Ketone readings deliver plasma equivalent results and are displayed in millimoles of β-Ketone per liter of blood (mmol/L).
The β-Ketone test measures β-Hydroxybutyrate (β-OHB), the most important of the three β-Ketone bodies in the blood. Normally, levels of β-OHB are expected to be less than 0.6 mmol/L.
β-OHB levels may increase if a person fasts, exercises vigorously or has diabetes and becomes ill. If your β-Ketone result is 0.0 mmol/L, repeat the β-Ketone test with new test strips. If the same message appears again or the result does not reflect how you feel, contact your healthcare professional. Follow your healthcare professional's advice before you make any changes to your diabetes medication programme. If your β-Ketone result is between 0.6 and 1.5 mmol/L, this may indicate development of a problem that could require medical assistance. Follow your healthcare professional's instructions. If your β-Ketone result is higher than 1.5 mmol/L, contact your healthcare professional promptly for assistance. You may be at risk of developing diabetic ketoacidosis (DKA).[1-6]
Please consult your doctor to determine a target range that works best for you.

Questionable or inconsistent results
If your test results are unusual or inconsistent with how you are feeling:
► Make sure the confirmation window of the test strip is completely filled with blood. ► Check the expiry date of the test strips. ► Check the performance of your meter and test strip with the control solutions. ► Make sure your monitor has correct coding, and the code is the same as the code.
Please Note: Unusually high or low β-Ketone levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

Quality Control Testing

Our control solutions contain a known amount of β-Ketone that can react with test strips. You can check the performance of the meter, test strip and your technique by comparing the control solution results with the range printed on the label of vial or individual foil pack. Checking regularly can ensure your test results are accurate. Please refer to the Owner's Manual for complete testing instructions.

IMPORTANT: The reference range of the control solutions may vary with each new test strip. Make sure you check the range on the label of vial or individual foil pack of your current test strip.

Chemical Components

> β-Hydroxybutyrate Dehydrogenase (*Pseudomonas sp.*) ≥ 0.5 U
> Mediator 55%
> NAD ≥ 0.5 µg
> Enzyme protector 8%
> Non-reactive ingredients 29%

Additional Information for Healthcare Professionals

Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples. Use fresh whole blood samples only. Professionals may use test strips to test capillary and venous whole blood.
Sample Size: 0.8 µL
Reaction Time: 10 seconds
System Measurement Range: 0.1 mmol/L to 8 mmol/L
Hematocrit Range: 10 % to 70 %

Accuracy
The table below displays how often FORA achieves this target. The chart is based on a study carried out on 160 patients (each patient was tested three times which resulted in 480 test results) to see how well FORA performed compared to Randox D-3 Hydroxybutyrate (Ranbut) reagent reference method results.

Capillary samples (n=480)	β-ketone concentration (mmol/L)	Regression analysis
	Range: 0.10 to 6.75 Mean: 1.13	y = 0.9245x - 0.0057, R^2 = 0.9335
Venous samples (n=480)	β-ketone concentration (mmol/L)	Regression analysis
	Range: 0.10 to 6.75 Mean: 1.12	y = 0.9926x - 0.0554, R^2 = 0.9772

Precision
In repeatability test, the standard deviation (SD) is within 0.1 mmol/L for each β-Ketone concentration < 1 mmol/L and the coefficient of variation (CV) is less than 7.5% for each β-Ketone concentration ≥ 1 mmol/L.
Use only with FORA 6 Duo / FORA 6 Plus / FORA 6 Connect Multi-functional meter.

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Table A summarized the substances tested at recommended concentrations which bias is less than 10% and demonstrated no interference with FORA Test Strip.

Substance	Concentration Tested
Captopril	23 µmol/L
L-DOPA	50 mg/L
Dopamine	5.87 µmol/L
Gentisic acid	117 µmol/L
Paracetamol	1324 µmol/L
Uric acid	3 mmol/L
Ascorbic acid	228 µmol/L
Unconjugated bilirubin	400 µmol/L
Cholesterol	15 mmol/L
Triglycerides	30 mmol/L

Reference
[1] Schade DS, Eaton RP. Metabolic and clinical significance of ketosis. Special Topics in Endocrinology and Metabolism 1982;4:1-27. [2] Wiggam MI, O'Kane MJ, Harper R, Atkinson AB, Hadden DR, Trimble ER, Bell PM. Treatment of diabetic ketoacidosis using normalization of blood β-Hydroxybutyrate concentration as the endpoint of emergency management. Diabetes Care 1997;20:1347-52. [3] Harano Y, Kosugi K, Hyosu T, Suzuki M, Hidaka H, Kashiwagi A, Uno S, Shigeta Y. Ketone bodies as markers for Type 1 (insulin-dependent) diabetes and their value in the monitoring of diabetes control. Diabetologia 1984;26:343-8. [4] Ubukata E. Diurnal variation of blood β-Ketone bodies in insulin-dependent diabetes mellitus and noninsulindependent diabetes mellitus patients: The relationship to serum C-peptide immunoreactivity and free insulin. Ann Nutr Metab 1990;34:333-42. [5] Luzi L, Barrett EJ, Groop LC, Ferrannini E, DeFronzo RA. Metabolic effects of low-dose insulin therapy on glucose metabolism in diabetic ketoacidosis. Diabetes 1988;37:1470-77. [6] Hale PJ, Crase J, Nattrass M. Metabolic effects of bicarbonate in the treatment of diabetic ketoacidosis. Br Med J 1984;289;1035-8.

Symbol / Symbol / Symbolel / Símbolo / Symbool / Símbolo / Símbolo / Symbol	Referent / Bedeutung / Signification / Riferimento / Omschrijving / Referência / Referencia / Znaczenie
IVD	• *In vitro* diagnostic medical device • Für die *In-Vitro*-Diagnostik • Dispositif médical de diagnostic *in vitro* • Dispositivo medico per uso diagnostico *in vitro* • Dispositivo médico para diagnostico *in vitro* • De unică folosință • Dispositivo Médico para Diagnostico *in vitro* • Urządzenie medyczne do diagnostyki *In vitro*
	• Do not reuse • Nicht wiederverwenden • À usage unique • Non riutilizzare • Niet voor hergebruik • Não reutilizável • No reutilizable • Nie używać ponownie
	• Consult instructions for use • Gebrauchsanweisung beachten • Consulter le mode d'emploi • Leggere le istruzioni per l'uso • Raadpleeg instructies voor gebruik • Consultar instruções de utilização • Consultar instrucciones de uso • Sprawdź instrukcję obsługi
	• Temperature limitation • Temperaturbegrenzung • Plage de températures limites • Limite di temperatura • Temperatuurbegrenzing • Limitação de temperatura • Limitación de temperatura • Ograniczenie temperatury
	• Use by • Zu verwenden bis • Date limite d'utilisation • Data di scadenza • Houdbaarheidsdatum • Utilizar antes de • Utilizar antes de • Data przydatności
REF	• Catalogue number • Katalognummer • Numéro de catalogue • Numero di catalogo • Catalogusnummer • Número de catálogo • Número de catálogo • Numer katalogowy
6M	• Use within 6 months after opening • Innerhalb von 6 Monaten nach dem Öffnen verwenden • Utiliser dans les 6 mois après l'ouverture • Usare entro 6 mesi dall'apertura • Te gebruiken binnen 6 maanden na opening • Utilizar no prazo de 6 meses após a primeira abertura • Una vez abierto, utilizar durante los primeros 6 meses • Po otwarciu zużyć w ciągu 6 miesięcy
CE 0123	• CE mark • CE-Zeichen • Marquage CE • Marcatura CE • CE-keurmerk • Marca CE • Marca CE • Znak CE
LOT	• Batch code • Chargencode • Code du lot • Codice di lotto • Partijcode • Codul lotului • Código de lote • Código de lote • Kod partii
	• Humidity limitation • Einschränkungen bei Luftfeuchtigkeit • Limite d'humidité • Limite di umidità • Vochtigheidsbeperking • Limitação de umidade • Limitación de humedad • Ograniczenie wilgotności
	• Do not use if package is damaged • Bei beschädigter Verpackung nicht verwenden • Ne pas utiliser si l'emballage est endommagé • Non utilizzare se la confezione è danneggiata • Niet gebruiken als de verpakking beschadigd is • Não utilizar se a embalagem estiver danificada • No utilizar si el envase está dañado • Nie używać, jeśli jest uszkodzone opakowanie
	• Manufacturer • Hersteller • Fabricant • Produttore • Fabrikant • Fabricante • Fabricante • Producent
#	• Model number • Modell-Nr. • Numéro de modèle • Numero modello • Modelnummer • Número do modelo • Número de modelo • Numer modelu
Σ	• Quantity • Menge • Quantité • Quantità • Aantal • Quantidade • Cantidad • Ilość
EC REP	• Authorized representative in the European Community • Autorisierter Bevollmächtigter für die EU • Représentant autorisé de la Communauté européenne • Rappresentante autorizzato per l'Unione Europea • Erkende vertegenwoordiger in de Europese Unie • Representante autorizado na União Europeia • Representante Europeo Autorizado • Autoryzowany przedstawiciel w Unii Europejskiej
UDI	• Unique device identification • Einmalige Produktkennung • Identifiant unique des dispositifs • Identificazione unica del dispositivo • Unieke apparaatherkenning • Identificador de dispositivo único • Identificador único del producto • Unikalny identyfikator urządzenia

For self-testing / Für Selbsttests / Pour auto-mesure / Per automisurazione / Voor zelf testen / Para auto-teste / Para automonitoreo / Do samodzielnego wykonywania badań

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EC REP **MedNet EC-REP GmbH**
ForaCare Suisse AG
Smart OTC GmbH