Official Website

JUMPER Health

Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan,Shenzhen, Guangdong, China

Tel: +86 755-26696279 **E-mail:** info@jumper-medical.com **Web:** www.jumper-medical.com

PRODUCT BROCHURE



SINCE

1998

JUMPER

Healthcare within mind

www.jumper-medical.com



TENS THERAPY DEVICE



APP Control · 5 Modes 15 Levels · One Button Control · Automatic Shutdown · Relieve Pain · Rechargeable Battery

RELIEVE PAIN, SOOTHE MIND.

Small, portable and powerful TENS devices from JUMPER offer drug-free pain relief for sore or aching muscles of the lower back, shoulder, arm, leg and more.

The unit also features one or two independent stimulation channels and four self-adhesive electrodes. It offers a wide range of functions for increasing general well-being, pain relief, maintaining physical fitness, relaxation, muscle revitalization and combating tiredness. For these purposes, you can either choose from pre-set programs or specify your own to suit your individual needs.

FEATURES

- Pain relief
- Pad off Alert
- Automatic shutdown
- Low battery indication
- Ultra-thin & lightweight
- Easy to charge and use
- Various massage methods
- Powerful & effective program





JPD-ES100

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
1-60 mins	1-120HZ	0-15

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
30V at 500Ω	LITHIUM	0.15-0.20ms

DIMENSION(mm)
MAIN UNIT 44 x 44 x 13
ELECTRODE PAD φ90 / φ60

MAIN UNIT WEIGHT
50g

APP Control · One Button Control · Pad Off Alert · 5 Modes · Rechargeable Battery · Low Battery Indication



A

JPD-ES210

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-120HZ	0-15

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.1-0.2ms

DIMENSION(mm)
MAIN UNIT 132.3 x 64.8 x 19.6
ELECTRODE PAD 98 x 64

MAIN UNIT WEIGHT
105g (WITHOUT BATTERY)

14 Modes · 15 Levels Intensity · Pad Off Alert · Available For 2 Users · Low Battery Indication · Automatic Shutdown



JPD-ES200

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-120HZ	0-15

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.15-0.20ms

DIMENSION(mm)
MAIN UNIT 121.8 x 63.6 x 22.2
ELECTRODE PAD 50 x 50

MAIN UNIT WEIGHT
86g (WITHOUT BATTERY)

5 Modes · 15 Levels Intensity · Pad Off Alert · Available For 2 Users · Low Battery Indication

NEW ARRIVAL

TENS
MASSANT
EMS



B

JPD-ES210

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-15

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
48mA at 500Ω	DRY BATTERIES	0.05-0.30ms

DIMENSION(mm)
MAIN UNIT 132.3 x 64.8 x 19.6
ELECTRODE PAD 98 x 64

MAIN UNIT WEIGHT
105g (WITHOUT BATTERY)

41 Programs · 15 Levels Intensity · Pad Off Alert · Available For 2 Users · Low Battery Indication · Automatic Shutdown

NEW ARRIVAL



JPD-ES220

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	1-100HZ ± 10%	0-20

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	80μs

DIMENSION(mm)
MAIN UNIT 138x 68.7 x 13
ELECTRODE PAD 50 x 50

MAIN UNIT WEIGHT
123.9g

Type-C · Rechargeable Battery · 20 Modes · 20 Levels Intensity · Available For 2 Users · Pad Off Alert

NEW ARRIVAL

TENS
MASSAGE
EMS



JPD-ES230

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-40

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	DRY BATTERIES	0.05-0.30ms

DIMENSION(mm)
MAIN UNIT 130 x 64 x 23
ELECTRODE PAD 50 x 50

MAIN UNIT WEIGHT
98g (WITHOUT BATTERY)

44 Programs · 40 Levels Intensity · Pad Off Alert · Available For 2 Users · Low Battery Indication · Automatic Shutdown

NEW ARRIVAL

TENS
MASSAGE
EMS



LPA11

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-100HZ	0-40

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	DRY BATTERIES	0.05-0.30ms

DIMENSION(mm)

MAIN UNIT 138x 68.7 x 13

ELECTRODE PAD 50 x 50

MAIN UNIT WEIGHT
123.9g

70 Programs · 40 Levels Intensity · Pad Off Alert · Available For 2 Users · Low Battery Indication · Automatic Shutdown

NEW ARRIVAL

TENS
MASSAGE
EMS



LPA10

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-40

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	0.05-0.30ms

DIMENSION(mm)

MAIN UNIT 151 x 78 x 32

ELECTRODE PAD 50 x 50

MAIN UNIT WEIGHT
145g
(WITHOUT BATTERY)

44 Programs · 40 Levels Intensity · Type-C · Pad Off Alert · Available For 2 Users · Low Battery Indication

NEW ARRIVAL

TENS
MASSAGE
EMS



LPB10

TREATMENT TIME	PULSE FREQUENCY	INTENSITY
5-60 mins	2-120HZ	0-40

MAX OUTPUT INTENSITY	POWER SUPPLY TYPE	PULSE WIDTH
50V at 500Ω	LITHIUM	0.05-0.30ms

DIMENSION(mm)

MAIN UNIT 132.3 x 64.8 x 19.6

ELECTRODE PAD 50 x 50

MAIN UNIT WEIGHT
105g
(WITHOUT BATTERY)

44 Programs · 40 Levels Intensity · Type-C · Pad Off Alert · Available For 2 Users · Low Battery Indication

TENS Therapy Device Parameters		 A		
Product Model	JPD-ES220	JPD-ES210	JPD-ES200	JPD-ES100
TENS EMS MASSAGE	TENS			
Treatment Time	5-60 mins	5-60 mins	5-60 mins	1-60 mins
Power Supply Type	Lithium	Dry batteries	Dry batteries	Lithium
Max Output Intensity (at 500Ω)	50 V	48 mA	48 mA	30 V
Pulse Width	80μs	0.1-0.2ms	0.15-0.20ms	
Pulse Frequency	1-100Hz	1-120Hz		
Low Battery Indication	√			
Number of Channels	2			1
Auto Power Off	√			
Bluetooth	—			✱
Main Unit Size (mm)	138 x 68.7 x 13	132.3 x 64.8 x 19.6	121.8 x 63.6 x 22.2	44 x 44 x 13
Electrode Pad Size (mm)	50 x 50	98 x 64	50 x 50	φ90 / φ60
Main Unit Weight (g)	123.9	105	86	50
Modes	20	14	5	

3 in 1 Combo Electrotherapy Device

Parameters



Product Model	JPD-ES210	JPD-ES230	LPA11	LPA10	LPB10
TENS EMS MASSAGE	TENS - EMS - MASSAGE				
Treatment Time	5-60 mins				
Power Supply Type	Dry batteries			Lithium	
Max Output Intensity (at 500Ω)	48 mA	50 V			
Pulse Width	0.05-0.3ms				
Pulse Frequency	2-120Hz				
Low Battery Indication	✓				
Number of channels	2				
Auto Power Off	✓				
Bluetooth	—				
Main Unit Size (mm)	132.3 x 64.8 x 19.6	130 x 64 x 23	113 x 61 x 22	151 x 78 x 32	151 x 82 x 33
Electrode Pad Size (mm)	98 x 64	50 x 50	50 x 50	50 x 50	50 x 50
Main Unit Weight (g)	105	98	95	145	128
Programs	41	44	70	44	44

