BESTMAN

ISO 13485 CE ISO 9001

BF-500++





BESTMAN

Ultrasound Fetal Dopple

BESTMAN

BF-500++

20 years experience of manufacture fetal Doppler, Quality Promised



individual light probe

gentle outlook

simaultanous Detection

product comparison

2.5 mhz probe



Product brand	BESTMAN	Battery	Built-in 9.6 v
Ultrasonic Frequency	2.5MHz±10%	Isptp	⩽0.1MPa
Atmospheric pressure	86kPa～106kPa	Ultrasonic intensity	⩽10mW/cm²
FHR range	60bpm～210bpm	Comprehensive Sensitivity	⩾90dB
Adapter input voltage	AC 220v±10% 50Hz±1Hz	Weight	⩽0.1MPa

Charger	X1
Hanging belt	X1
Manual / warranty card	X1
certificate of quality	X1

Security certification

BESTMAN fetal doppler, a well-known brand

ISO 13485 CE ISO 9001

BESTMAN

BF-600+

OPERATION GUIDE

Pick up the machine and gently push the probe back to remove the probe



After removing the probe and apply 50g coupling agent to the probe



According to the time of pregnancy to detect fetal heart





CORRECT USE OF SKILLS



Proper way of holding



Correct use of angle

SHENZHEN BESTMAN INSTRUMENT CO.,LTD

www.szbestman.cn

Add: 8/F,YiFang Bldg,No.315 Shuangming Ave,Dongzhou Community, Guangming St.,Guangming Dist., Shenzhen, Guangdong,518107, China

Tel: +86-755-26411423

Fax: +86-755-26495167