

PROFESIONAL URBAN FIRST AID

SAVE LIFE IN 4 MINUTES

GERMAN TECHNOLOGY

ONEMED Automatic External Defibrilator


Oneheart B8

ONEMED


TO SAVE LIVES EVERYWHERE


Cardiac arrest can happen to anyone, anywhere, anytime

Sudden cardiac death (SCD) refers to the sudden stop of the pumping function of the heart in any patient at an unexpected period of time, resulting in interruption of blood flow to the brain and various organs of the whole body, loss of consciousness, respiratory arrest, and even sudden death.



THE ONLY FIRST AID MEASURES IN CARDIAC ARREST: CPR AND DEFIBRILLATION





CONVENIENT AND EASY TO USE

Everyone can save lives if everyone can use first aid equipment







•Short operation time:
turn on immediately by opening the cover, the operation interface is focused with visualization, less thinking time required.

•One-key switch:
with one-key switch between adult/child mode and language mode, the operation becomes more convenient.

•Double reminder:
double reminder with button sound control and light control, the operation becomes safer.

•Full intelligent guidance for CPR:
The CPR feedback device can provide real-time feedback on the quality of CPR, and real-time encouragement and feedback of the CPR situation through voices such as "press hard" and cycle sounds, making first aid simpler


EXTREMELY FAST AND EFFICIENT


In a race against time, technology makes first aid more efficient

The Oneheart B8 AED adopts pre-charging technology, and the time from power-on to discharge-ready is only 6 seconds.



The Oneheart B8 AED adopts advanced intelligent biphasic waveform defibrillation technology, which has a high success rate in first try of defibrillation and with little myocardial damage. The corresponding energy can be automatically configured according to the patient's own impedance, which is suitable for a wider group.






RELIABLE AND DURABLE

Adapt to various harsh environments and can be used in a variety of scenarios

The Oneheart B8 AEDs have passed more than 200 quality tests before leaving the factory. After leaving the factory, self-inspection is carried out daily, weekly, monthly and semi-annually, including the expiration date and status of each module of the device, and the electrode pads and the batteries. Once a fault occurs, the status bar can display a warning.




• The working temperature of the series AEDs is -5°C~+55 ℃/23 ℉~131 ℉. After entering a --20°C environment from one at room temperature, the device can work for at least 20 minutes, and the device can be applied in various harsh environments.



• The series AEDs have a protection level of IP55 and a drop-resistant height of up to 1.5m, and are suitable for ambulance transport. The whole device has no exposed interface (except the electrode pad interface), and can be used in various scenarios.



• Data can be obtained through USB mobile disk, which is convenient, fast, safe and reliable.



Simple operation | Smart device | Data transfer | Voice broadcast | Fault warning | Equipment self-check





model	Oneheart B8
language (voice)	4
CPR feedback sensor	●
Status screen	●
CPR press prompt	●
From AED analysis to charge 200J done	≤5s
From AED analysis to charge 360J done	≤12s
From open device cover to charge 200J done	≤7s
From open device cover to charge 360J done	≤14s
Energy output	Default adult energy sequence : Level 1:200J Level 2:300J Level 3:360J Default children energy equence : Level 1:50J Level 2:70J Level 3:100J
Data storage	ECG wave: 10h Event: 2000 Audio log: 2h CPR data: 10h Self-test report: 3650
Data transfer	USB2.0
Buttons	Power、Child、Language
Battery type	LiMnO2 ,24V, 1.2Ah, non-rechargeable
Pad type	Disposable defibrillation electrodes foradult and pediatric, default without CPR feedback sensor
Dimensions	22 x 9.7 x 29.6 cm (L x W x H)
Weight	approx. 2.5kg (with energy module and pads)
Shelf life	10years

JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD
Manufacturing Address: No.1 Baisheng Road, Development Zone, Danyang, Jiangsu Province, P.R.China
Tel: 86-0511-8690 0859 86-0511-8690 0860
E-mail: yuwellsales@yuyue.com.cn; sales@yuyue.com.cn
www.yuwell.com

KEMENKES RI AKL 30505420031
PT. INTISUMBER HASIL SEMPURNA GLOBAL
Surabaya, Jawa Timur, Indonesia
www.onemed.co.id