stomdevice

BESTMAN®

深圳市贝斯曼精密仪器有限公司
SHENZHEN BESTMAN INSTRUMENT CO.,LTD

关注健康 共赢未来





The destruction temperature of the hypodermic needle is 1200 degrees.

Syringe Destroyer BD-310

Feature:

Melt the needle and cut the syringe head and infusion tube, alarm prompts and ozone disinfection;

Destroy needles size rangement:18G-34.

Specification:

Power supply: AC 220V±10% 50Hz±lHz AC 110V±10% 60Hz±lHz

Input power: BD-310<80W

Working temperature : +10°C~+40°C Fuse specification:

(AC! 220V) BD-310 250V/2 AL (ACI 110V) BD-310 125V/4 A L

Quantity melted per time: >10000 syringes

The destruction temperature of the hypodermic needle is 1200 degrees.

N.W. (KG): 2.36

Dimension(CM):19.6x23x14

Add: 8/F, Yi Fang Building, No. 315 Shuang Ming Blvd, Guangming Dist., Shenzhen, P.R. China (Factory)
Tel: +86-755-2671 3783; 2671 3784 Fax: +86-2649 5167; 2671 3783; 2671 3784